Protocole :

	Ø	%	L	Stop	Bague	RPM	Recommandation Couple (N.cm)
Coronal	25	10	17	-	○	300 - 400	2 - 3
Median	25	6	21-25	○	●	300 - 400	2 - 3
Apical1	20	4	21-25-29	○	●	300 - 400	2 - 3
Apical2	20	6	21-25-29	○	●	300 - 400	2 - 3
Apical3	30	6	25	○	●	300 - 400	2 - 3
Apical4	40	7	25	○	●	300 - 400	2 - 3
PreSequence	17	2	21-25	○	●	300 - 400	2 - 3

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
1. Indication

CMA: root canal preparation during non-surgical endodontic treatment and retreatment.

For use by dental professionals only.

2. Contraindications

There is no contraindication to the use of CMA products after the dental professional has decided to perform a non-surgical root canal treatment.

3. Complications

In cases of complex canal anatomy, per-operative risks (instrumental breakage, ledge, stripping, zipping, false path, perforation, etc.) could occur and lead to a risk of infectious processes.

4. Warnings/precautions

- The decision to use an endodontics instrument must be based on the clinical case expertise, particularly where the canal anatomy is considered to be too complex.
- Patients identified as having a risk of infectious endocarditis.
- Contains Nickel and Titanium and should not be used on patients with known allergic sensitivity to these metals.
- Respect the good dental practice, in particular the use of a dental dam and gloves.
- Use in continuous rotation at the recommended speed.
- In cases of complex anatomy, the maximum torque may differ from the value recommended by the manufacturer.
- Use according to the recommended procedure (§7).
- Instruments supplied non-sterile: Follow the reprocessing instructions (§8) before use.
- Keep the UDI information which only appears on the labelling until the last use.
- Year of manufacture: see labelling.
- Do not use if doubt on the respect of storage conditions.
- In case of doubt concerning the product identification, do not use.
- Check condition of the instrument before use between each canal. If the instrument is damaged or shows signs of wear, do not use it.
- Inform the manufacturer and the national regulatory authority of any serious incident relating to the instrument.

5. Clinical claims

CMA in normal use conditions.

Clinical performance:

- Shaping ability: adapted debris removal and cutting efficiency
- Respect of canal anatomy
- Geometry adapted to irrigation*
- Maintains the apical foramen in its original position**

Safety :

- Limited risk of breakage
- Limited peroperative risks*
- Limited risk of screwing effect

*Not applicable for Coronal and PreSequence.

**Not applicable for Coronal.

6. Characteristics

Type 2 rotating instruments under EN ISO 3630-1: 2019 standard:

- Maximum recommended torque (not required): see table 'Recommended Torque'.

Material of operative part: Nitinol and stainless steel (for PreSequence instrument).

Use in connection with an endodontic contra-angle according to EN ISO 1797: 2017 (Type 1).

Number of uses:

- maximum of 10 canals recommended if the instrument is not visibly damaged.

- maximum of 3 canals recommended for PreSequence file.

Medical device class according to Directive 93/42/EEC and MDR 2017/745: Ila.

7. Protocol

See Protocol.

8. Instructions for use

General recommendations	- For all metal devices, the use of anticorrosion disinfecting and cleaning agents is recommended. - For your own safety, please wear personal protective equipment (gloves, glasses and mask). - Do not use cleaning or disinfecting agents containing phenol, aldehyde and alkaline composition. - Always respect the instructions for use provided by the manufacturer of the products.
-------------------------	---

Reuse limitations	- Due to the product design and the materials used, no definite limit to the maximum number of performable processing cycles can be specified. The service life of the medical devices is determined by their function and careful handling. Multiple use, disinfection and re-sterilisation cycles may lead to increased risk of file separation. - The user must ensure that the processing method used, including resources, materials and personnel, is appropriate and meets the applicable requirements. - The state of the art and national laws require that validated procedures be followed.
Material needed	- Gloves, masks, gown as recommended by the manufacturer of cleaning agent and detergent - Tap or deionized water* - Disinfectant (neodisher® Septo Active) - Detergent (neodisher® MediZym) - Small soft brushes - Container - Ultrasonic tub or washer-disinfector - Class B sterilisation apparatus <i>Remark: All material used should be cleaned and replaced regularly. Identify material used for each step of the process (initial treatment, cleaning or rinsing).</i>

1	Initial treatment Place used products in a container or a wipe with tap water at 20-40 °C and 1.0 % neodisher® Septo Active for 5 to 15 min. Rinse the products with tap water at 20-40 °C for 1 min. If there is waiting time before the next step, make sure the device stays moist by placing it in a wet wipe. Do not exceed 1 hour waiting time. <i>Remarks:</i> - Do not use fixing agents or hot water (> 40 °C), as this causes fixation of residues and can impair successful cleaning. - Follow instructions and respect concentrations and immersion times given by the manufacturer (an excessive concentration may cause corrosion or other defects on devices).
2	Preparation before cleaning If the devices have visible impurities, it is recommended to manually pre-clean by brushing under tap water at 20-40°C for at least 1 min with a soft brush until all impurities have been removed. <i>Remark: Follow instructions and respect concentrations and immersion times given by the manufacturer (an excessive concentration may cause corrosion or others defects on devices).</i>
3	Visual inspection Inspect the used products and discard damaged products (broken, untwisted or abnormally bent).
4	Cleaning Insert the products into an ultrasound apparatus beaker. Run ultrasound apparatus for 10-30 min with tap water and 0.5-2.0% neodisher® MediZym. <i>Remarks: Follow instructions, observe water quality, concentrations and cleaning time specified by the manufacturer of the cleaning solution.</i>
5	Rinsing Rinse the products with tap water at 20-40°C for 1 min. <i>Remark: It is recommended to use deionized water.</i>
6	Drying Dry the products with compressed air until products are visibly dry.

automatic	4 5 6	Cleaning/Rinsing/Drying Place the instruments in the tray of the washer/ disinfectant's sliding trolley. Perform cleaning cycle with 0.2-1.0% neodisher® MediZym. Perform drying. <i>Remarks:</i> - Disinfection (thermal or chemical-thermal) is not required since the products are sterilised after cleaning. - Follow instructions and concentrations specified by the manufacturer of the detergent solution. - Follow the instructions of the washer-disinfectant and verify the success criteria have been met after each cycle as specified by the manufacturer. - The final rinse step should be with deionized water. For other steps, follow the water quality defined by the manufacturer. - Use only approved washers-disinfectors that comply with EN ISO 15883, and are maintained and validated regularly.
7		Visual inspection Inspect the used products. Re-do steps 4,5 and 6 if the product is visibly not clean or discard any damaged products.
8		Packaging Place the instruments in a paper-plastic pouch for steam sterilisation in compliance with ISO 11607 and EN 868 standards. <i>Remarks:</i> - For sharp devices that are not contained within a box, silicon tubes should be placed around the devices to prevent the packaging from being pierced. - Seal the pouches as recommended by of the pouch manufacturer. If a thermo-sealer is used, the process must be validated and the thermo-sealer must be calibrated and qualified.
9		Sterilisation In accordance with French regulations: - Apparatus: class B - Minimum temperature: 134°C - Minimum time: 18 min - Absolute pressure: 2.2 bar - Minimum drying: 20 min Other regulations: - Apparatus: class B - Minimum temperature: 132°C - Minimum time: 3 min - Absolute pressure: 2.2 bar - Minimum drying: 20 min <i>Remarks:</i> - When sterilising multiple instruments in one autoclave cycle, ensure that the steriliser's maximum load is not exceeded. - Place the pouches in the steam steriliser according to the recommendations given by the steriliser manufacturer. - Use only a pre-vacuum air removal steam steriliser that meets the requirements of EN 13060 (class B, small steriliser) and EN 285 (full-size steriliser), with saturated steam.
10		Storage after sterilisation Store the products in a dry, clean and dust-free environment at the temperature specified by the paper-plastic pouch by the steam steriliser manufacturer. <i>Remark: Check the packaging and the medical devices before using them (packaging integrity, no humidity and use-by date). In case of damage, a complete rerun should be performed.</i>

9. Storage and transport conditions

Keep products away from heat and moisture.

10. Disposal

After use, instruments must be placed in a secure container, used to collect cutting or sticking instruments (like needles or disposable bistouries) as per good dentistry practice.

11. Symbols


	NiTi raw material		Quantity
	Stainless steel material		Assortment
	'CLASSICS' handle		Do not use if packaging is damaged
	Manufacturer		Non-sterile
	Batch code number		Date of manufacture
	Sterilizable in a steam sterilizer (autoclave) at the temperature specified		
	Caution		Medical device
	Keep dry		Keep out of direct sunlight

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